

# CASE STUDY

**LEGACY CLEANING PROCESS OPTIMIZATION** 

# LEGACY IMPROVEMENT PROCESS IMPROVEMENTS HYDE CORE LAB

Dublin, Ireland – MSD Swords Manufacturing Facility
United States – Hyde CORE Laboratory – Louisville, Colorado

#### PROJECT DESCRIPTION

A legacy cleaning process for the removal of a product residue containing significant amounts of the Active Pharmaceutical Ingredients (API) residue was identified by the engineering team as both a potential safety hazard and a known bottleneck in the manufacturing process.

The cleaning process could take up  $_{to}$  five (5) working days for Cleaning in Place (CIP) operations. The long duration of the cleaning process was due to an inconsistent cleaning cycle that did not completely remove the hormonally bases oral solid dose residue from the process equipment. Due to the inconsistent nature of the residue removal, the process routinely required manual intervention, including confined space entry resulting in a breach of the room integrity. As the product residue represents the highest danger at a Level-5 on the Occupational Exposure Band (OEB), manual cleaning was not desirable.

## SCOPE OF SERVICES PROVIDED

### Mechanical

- Replaced existing spray balls with new Toftjorg spray ball design (Greater cleaning efficiency)
- New pipe work & spray ball locations in Bolz vessel based on computer imaging (Complete spray coverage)
- New 6kW Air Heater (Increased thermal mass)
- New Air supply header (Continuity of air supply)
- New Conductivity probe (Guarantee homogeneity of CIP solution)
- New detergent pumps & set-up (Obviates air locks & priming issues)

## **Automation**

- A new recipe designed from 1st principles to clean and dry
- New intelligent detergent dosing program incorporated

### Scale up of Development & Pilot CIP process

- New designed recipe, developed at Hyde CORE lab introduced to production scale
- Successful Riboflavin & "Edge of Failure" performed
- Installation, Operation, & Performance Qualification (IQ/OQ/PQ) as well as Commissioning Validation (CV) successfully complete on schedule/budget



Pharmaceutical ENGINEERING

Nov/Dec 2014, V 24, No 6, p50-60

"Application of Lean Six Sigma to Optimize a Legacy Cleaning Process."

By Emmet Manning, Brian Earls, Keith Bader, Fearghal Downey PhD,
and Kelly Scalva